

# How to Import Previous Database (M8000) Using MyoVision software Atlas & Apricot

Last Modified on 02/19/2018 4:35 pm PST

## Summary

If you recently upgraded from the M8000 software and are attempting to import a previous database file to your new installation of MyoVision Atlas or Apricot, you may encounter a Jet database engine error. Your database is saved in a database governed by Microsoft Access, which will hiccup and cause this error. The MyoVision software no longer uses this database format because of its tendency to cause errors. If you encounter this error you are using an older version of the software. This document will go into the specifics of resolving this error.

## Before You Begin


- You will want to make sure every step of the way that you have created a backup of your existing database.
- Backup your files before proceeding with any database recovery. There is no guarantee that your database will be 100% intact.

## 1. Open MyoVision Software

Open your **new installation of the Myovision Atlas or Apricot software**. Click on the icon your desktop or from your **Start** menu


## 2. Previous Database

Click **File** and select **Import** then click on **Previous Database**



### 3. Important Process

Select OK when you receive the Important Process warning prompt




### 4. Locate Database

You will be prompted to locate and select the database file you wish to


import. Please note that the default directory is C:/WinScan98. The file for the old Access database is **database.mdb**. *If your database is located in an alternate directory, you will need to locate where your database is stored in order to obtain the correct directory path.*

*Note: If two .mdb files are listed, take the larger size of the ones listed. This is your database; you want the most data available.*




## 5. Open Database

Click **Ope** to begin importing your previous database. Once you select **Open** the MyoVision software will begin analyzing your existing database.



## 6. Upgrade and Import

Once your existing database has been analyzed, the upgrade and import will begin automatically. Please note that depending on the size of your database, this process can take up to 1 hour.



## 7. Verify

Once the import process has completed, please select the **Open Exam** command from the MyoVision software to verify that your previous database has imported successfully.



## Getting Started

Select one of the options below or, above.

Start an Exam

Open

Select Patient Exam

Patient:			Exam:		
Name	ID:	Gender	Name	Date	Time
Be		M	Full Spine	Mar 31, 2001	02:37:43 PM
Be		F	Full Spine	Apr 07, 2001	04:12:12 PM
Be		M			
Be		M			
Be		F			
Be		F			
Be		F			
Be		M			
Be		M			
Be		M			
Be		F			
Be		M			
Be		M			
Be		M			
Be		M			
Be		M			
Be		M			
Be		M			
Be		M			
Be		F			
Be		F			
Be		M			
Be		F			
Be		F			

Database Path Source:  
Provider=SQLLEDB.1;Integrated Security=SSPI;Persist Security Info=False;Initial Catalog=MyoVision;Data Source=(local)\SQLEXPRESS;Connect Timeout=15

Comparison: Hold down <Ctrl> key, click Exams to compare, click Show Exam.

OK  
Cancel