How to Edit Patient Information

Last Modified on 07/09/2018 4:18 pm PDT

This article will show you how to edit a patient's name, ID, and gender.

1. Open the Patient menu and select Edit



2. Select the patient you want to edit and click Edit.

Name	ID:	Gender	Most Recent Exam Date		
Ideal, Results	1212	M	Nov 22, 2014		
panpic, JCarr	1212	M	Jan 22, 2016		
				•	
				2 E dit	ì

3. Enter in the correct patient information. Then, click OK.

New	Name
'atient Identificatior 0000000	Number:

4. Click YES to confirm the changes



You have successfully edited the patient's information.