

Using the Text to Mobile Feature

Last Modified on 09/24/2019 10:33 am PDT

Before you Begin

You must have entered and completed the Escan setup portion of the software setup before beginning the Text-to-Mobile.

[Click Here](#) for a guide on how to setup Escan.

In order to use this feature you will need Feature Pack G. You can download the most recent version of the MyoVision software [Here](#) .

To purchase Feature Pack G and to continue to get all of the new updates and features in the MyoVision software purchase your Annual Software Subscription on our [Store Site](#) today.

Summary

This article will guide you through using the Text to Mobile feature of the newest MyoVision software.


Some **Verizon customers** have reported having issues receiving an EScan Text-to-Mobile. This is due to a known issue with the Verizon server that manages their MMS gateway. They are aware of the issue and working on it. Some customers have reported having better success using Verizon's old MMS gateway, which is [@mypixmessages.com](mailto:mypixmessages.com) .

US-based customers can [Click Here](#) to be taken to a free tool that shows you the MMS Gateway for any US-based phone number.

[MMS Carriers and Gateways](#)

[Send Exam Via Text to](#)

[Add a Mobile](#)



International Carriers are handled by entering the Gateway for the carrier and then inputting the phone number.

Click on the buttons inside the tabbed menu for video guides on using the new Text to Mobile feature:

