

# SQL [DNETLIB] Error

Last Modified on 06/15/2018 11:55 am PDT

## Table of Contents

- [Summary](#)
- [Is SQL Server Installed?](#)
- [Is SQL Service Running?](#)

---

## Summary

The purpose of this article is to discuss how to resolve the "SQL [DNETLIB] [ConnectionOpen...]" error. This error occurs when the software doesn't have access to the service the SQL Database it is being pointed at or the service isn't running.



---

## Is SQL Server Installed?

MyoVision software uses SQL Server to save the patient database data.

Step 1: Check that SQL Server is  
Installed

Step 2: Install SQL Server  
2014

Step 3: Test  
Database

---

# Is SQL Service Running?

When troubleshooting an SQL [DNETLIB] error make sure to verify that SQL Server service is running.

Step 1: Open Windows Services	Step 2: Find and Check the SQL Service
Step 3: Check the Service Startup Type	Tutorial Video
	Network Databases

---