Troubleshooting Issues With Interpretation Reports

Last Modified on 10/30/2025 2:07 pm PDT

Table of Contents

- Summary
- Check for Updates
- Computer Resources
- Pre-Open Word
- Interpretation Report Templates Not Showing Up
- Other Troubleshooting Steps
- Troubleshooting the Online Interpretation Service

Summary

Return to Table of Contents

This article will take you through the steps of troubleshooting errors and issues when trying to create an Interpretation Report

Check for Updates

Return to Table of Contents

If there are pending Windows updates, it can cause an error when trying to make interpretation reports. To check if you are fully updated, go to the Microsoft support article here: https://support.microsoft.com/en-us/windows/install-windows-updates-3c5ae7fc-9fb6-9af1-1984-b5e0412c556a

You will want to make sure that the online service is enabled.

1. Locate the "toolbar" at the top of the software and left-click "Setup", go down to "Settings", and select "General".




2. Make sure that the check box next to "Check for and use internet services if available" has a check

mark. If it does not, left-click the box to place a check mark. Click, OK.




Now try sending an exam for online interpretation. Make sure that the exam that you want to send off for online interpretation is open.

1. Right-click the exam, move the mouse so that the cursor is over "Create Interpretation Report" and select "Online Service".



2. In the window that opens select whether you have "Already Purchased" your online interpretation service or whether you would like to "Purchase Now" your online interpretation service.



3. Fill-out the "Exam Interpretation Information" fields. When, you are ready, click "Send".

