SoloVision Not Powering On

Last Modified on 07/20/2018 2:43 pm PDT

<u>Summary</u>

If your SoloVision isn't powering on, follow the troubleshooting steps below

Check USB Cable and Port

The SoloVision receives power directly from the USB cable that connects to the PC. Try unplugging the SoloVision and plugging it into a different USB port.

If the chord has a kink in it, has wiring exposed, or is over 2 years old you will want to try replacing the USB cable. Click the link below to purchase one directly from the MyoVision Store.

<u>http://store.myovision.com/product-p/pb-200.htm</u>

If you are still experiencing issues call our support desk at **1-800-969-6961 Option 2** for support for additional troubleshooting.