

How Do I Install the MyoVision Software?

Last Modified on 10/28/2019 10:50 am PDT

Summary

This document shows you how to download and install MyoVision software for the first time.

Would you like personalize help or one-on-one training with your software? Visit our store website [Here](#) to purchase.

Before you Begin

1. Make sure the Windows account you are logged into has Administrator privileges
 - Check [Administrator Privileges](#) if you are unsure.
2. **IMPORTANT:** Make sure Windows 10 is fully updated before updating MyoVision.
 - To check if you are fully updated, go to the Microsoft support article here:
 - [Windows 10](#)
 - [Windows 11](#)
3. You will need at least 10 GB of free hard drive space on your computer to install or upgrade the MyoVision software.
4. Installation can take up to 45 minutes to complete. Do not restart your computer or disconnect the internet connection.

Information for New Users

If this is your first time downloading and installing the MyoVision software, you automatically get a 30-day free trial of the full MyoVision software. After 30 days, if you have not entered a Registration Key, the software will only allow you to access the basic scan features.

A new software registration key is included with a one year software subscription purchase, by purchasing a new system, or by re-certifying a used system. For more information on Registration Keys [Click Here](#).

Note: If you bought a system from MyoVision a one year software subscription is included. After that updates for bug fixes are free but to get the latest feature packs you will need to re-subscribe to the subscription. Click the link to purchase the [Software Subscription](#) plan.

Customers that Purchased a Used or Second-Hand Device

The MyoVision Wirefree System is a Class II Medical Device. To ensure the device meets "factory new" specifications, as required by the FDA, a re-certification process is required prior to our supporting the product. This will also ensure that you are notified of any legal updates and information that may arise.

If you purchased your MyoVision system second-hand and are having issues getting it set-up you will need to re-certify your devices before our Support Team is able to assist you. [Click Here](#) to re-certify your device. Customers that choose to purchase second-hand devices are responsible for contacting the seller for any necessary software registration information.

How to Download and Install the MyoVision Software

1. Download MyoVision Installer

- Download the software at: <https://www.myovision.com/download.html>
 - Click the "Download Now" button and **save** the file.
 - Internet Explorer or Microsoft Edge: When prompted, choose to **Save** the installer program.
 - Chrome: It will automatically **Save** the file and the progress will be shown at the bottom of the browser window.
 - FireFox: When prompted, choose to **Save the file**.

2. Open the Installer

- Open an Internet Browser Window.
 - Internet Explorer or Microsoft Edge: Click **Run** at the bottom of the browser window.
 - Chrome: Click the Installer Icon at the bottom of the browser window.
 - FireFox: Click on the Download Icon next to the bookmarks Icon at the top of the browser window, then click on the Installer File.
- Open via File Explorer.
 - Most browsers will save files to the downloads folder on your computer. Double click on the installer file to launch.

3. Install MyoVision

1. When you first open the MyoVision installer, click **Next** to proceed with the install.
2. For new MyoVision installations, you may need to install a few Program Prerequisites first. If so, this pop up; Click **OK**.
3. The installer will tell you which programs need to be installed. Click **Next**.
 - **Note:** Programs listed will vary depending on what is already installed on the computer.
4. Installation generally takes several minutes to complete, sometimes up to one hour.
5. Once the Prerequisite Programs are installed, the MyoVision installer will start. Click **Next** to continue.
6. Next, **choose** the option that best fits your install. Click on the first or second option, then click **Next**.
 - **Note:** You will not be asked for your Registration Key during the install so if you do not have it you can still choose the second option.
7. Select your location, then click **Next**.
8. The Install is ready to start. Click **Next** to proceed.
9. Once the MyoVision software is installed, a confirmation window will appear with the Database connection status displaying **Success**.
 - If you wish to open the software at this point, click the box next to "Run MyoVision Software now."
 - If your Database connection status does not say **Success**, run the MyoVision software and

check to see if an SQL Server error appears.

- If you get an SQL Server error, view our support article [Common SQL Server Errors and Solutions](#) or call Technical Support at 1-800-969-6961 Opt 2.

4. Set up MyoVision software

- Once the software install is complete, you will need configure your software. Follow the links below for help with customization and setup:
 - [How Do I Enter my Software Registration Key?](#)
 - [Configuring Your MyoVision Software for Single or Dual Scanners](#)
 - [Setting Up EScan](#)
 - [Customize the EScan Message](#)
 - [How Do I Enter My Thermography Calibration Values?](#)
 - [How Do I Set Up My Wired Thermography Devices?](#)
 - [Configuring a Network Database in Your Office for MyoVision](#)
-