

# SQL 10 Error Troubleshooting

Last Modified on 03/10/2026 8:56 am PDT

## Summary

This document shows you how to troubleshoot an SQL 10 error

## Before you Begin


1. Make sure the Windows account you are logged into has Administrator privileges
  - Check [Administrator Privileges](#) if you are unsure.
2. **IMPORTANT:** Make sure Windows OS is fully updated before continuing.
  - To check if you are fully updated, go to the Microsoft support article here:  
<https://support.microsoft.com/en-us/windows/install-windows-updates-3c5ae7fc-9fb6-9af1-1984-b5e0412c556a>

## How to Troubleshoot a SQL 10 Error

These steps are listed in order of likeliness. Completing all of them is not required. Once the error is resolved, no further troubleshooting should be necessary.

### 1. Check the Database Service Status

- For network database setups, see our support article [Resolving SQL Errors On a Network Database](#)
- Go to **Setup > System > Database** and check the **Status** section circled below



- If the status shows:
  - **"Not Installed"**,
    - Most likely, SQL was not installed at the same time as MyoVision. **Make sure to update Windows** and follow the instructions in the [Installing SQL](#) support article.
    - If a database was previously on the computer (able to perform and review exams), see step 2.
  - **"Stopped"**,
    - If the MyoVision software is newly installed and has had no exams performed yet, see our support article [SQL Error code 10 - Service not running on New MyoVision software installation](#)
    - If the MyoVision software worked previously and had patient exams performed, perform the steps in our [SQL Error Code 10 - Service not Running](#)
      - If you already performed the steps in the SQL Error code 10 article linked above, move onto step 2. on this article

-  the software is configured for a network database, which means the database is on a different computer.
  - See our article [Resolving SQL Errors On a Network Database](#)
  - Not sure if you have a network database? Contact our technical support at 1-800-969-6961, option 2 (8:30am - 4:00pm Pacific Time, Monday - Friday) (+1 206-448-3464 International)

## 2. Check the Instance Name

- If you have multiple SQL instances installed, make sure the instance name is correct by following the instructions in the [Wrong SQL Instance](#) support article. If you have a network database setup, be sure to follow that set of instructions as the name will be different.

## 3. Check for installed AntiVirus software

- Popular antivirus software such as; Avast, AVG, McAfree, Norton, and Karpenski are all known to cause an SQL error code 10 or SQL error code timeout when their firewall systems are active.
  - Lookup and follow how to turn off the active protection and firewalls for your antivirus software.

**Note:** If you are still having trouble getting your MyoVision database to work properly, please feel free to contact our technical support at 1-800-969-6961, option 2 (8:30am - 4:00pm Pacific Time, Monday - Friday) (+1 206-448-3464 International)

---