

# How to Configure Software for Series-5 After Upgrading from a 3G/4G System

Last Modified on 12/18/2025 3:44 pm PST

## Summary

This document shows you how to configure your MyoVision software for use with a SERIES-5 MyoVision system.


## Before you Begin

1. Make sure your computer is updated to Windows 11. You can check if your computer has Windows 11 here: <https://www.microsoft.com/en-us/windows/get-windows-11?r=1>

## How to Configure MyoVision Software for SERIES-5

### 1. Use App Setup to Configure your software


- In the MyoVision Software, go to **Help > App Setup**. The following window should appear:
  - Click Next to begin setup and move through each selection window



- Select the VisionHUB option and **deselect the 3G/4G option** by clicking on it before clicking "Next".
  - **If one or more of your devices (other than Thermography) still use the 3G/4G hardware**, then leave the 3G/4G option selected and the software will prompt to choose the specific system hub for each device type.
    - Thermography device settings are stored separately and are not affected by the system hub setup



- Choose your devices accordingly.
  - If you are unsure whether you will be using 1 or 2 ScanVision, choose the option for "One Scanner". You can always go back and change the settings later. See our support article; [Configuring Your MyoVision Software for Single or Dual Scanners](#)



- When finished, click '**Done**'. Your software is now configured!