

Configuring Gmail to use with EScan

Last Modified on 03/28/2023 9:16 am PDT

Summary

This document shows how to setup your Gmail account with 2-step verification and an app password for use with EScan.

- What is 2-step verification?
 - When logging in to the Gmail account for the first time on a new computer, Gmail will send a text message to the attached phone number with a verification code that will need to be entered within 60 seconds. This is done to increase the security of the email account. After this, you can log in as normal.
- Is this really required? What if the email I want to use is used by several of my staff?
 - Yes, unfortunately Gmail announced mid-March 2022 that they are removing the secondary option (Less Secure App Access) on May 30st, 2022 and Gmail will block emails sent from MyoVision.
 - If you have a shared email address, we recommend creating a new email address for MyoVision to use. You can set the email to forward to your normal email address.
<https://support.google.com/mail/answer/10957?hl=en>


Before you Begin

1. Ensure your software is setup correctly to send emails. See our support article [Setting Up Escan](#)
2. If you have a Gsuite account contact MyoVision customer support for assistance at 1-800-969-6961 opt 2 (International: 206-448-3464) Office hours are 8:00am - 4:00pm M-F, PST.

How to Configure Gmail to use with EScan


1. Access Gmail Account Options

- Log in to the Gmail account being used. Click the 3x3 grid in the top right then click "Account"




2. Access Security Settings

- On the left-hand tab, click Security



3. Turn on 2-Step Verification

- If 2-Step Verification is already on, skip to Step 4
- Click on **2-Step Verification > Get Started**
- Follow the on-screen steps



4. Create an App-Password

- Under "Signing in to Google" select "2-Step Verification" again
- At the bottom, select "App passwords"

← 2-Step Verification



Security Key

A security key is a verification method that allows you to securely sign in. These can be built in to your phone, use Bluetooth, or plug directly into your computer's USB port.

Devices that don't need a second step

You can skip the second step on devices you trust, such as your own computer.



Devices you trust

Revoke trusted status from your devices that skip 2-Step Verification.

[REVOKE ALL](#)

App passwords

App Passwords aren't recommended and are unnecessary in most cases. To help keep your account secure, use "Sign in with Google" to connect apps to your Google Account.



App passwords

1 password

- Choose "Select App", click "Other" and call it "MyoVision" and click "Generate" and copy the generated 16 digit App Password
- Open MyoVision, go to **Setup > System > Email** and enter the generated App Password into the password field
 - **IMPORTANT:** If the Gmail password is reset or changed then any App Passwords will be deleted. A new App Password will need to be created and copied into MyoVision.

5. Send A Test Email

- Click "Send Test" and enter any email address that you have access to
- After pressing send it will show if the test succeeded
- If the email was unsuccessful, check to make sure that your password and email service provider is correct. If they are then contact MyoVision Technical Support for assistance
 - **IMPORTANT:** With Gmail accounts do not use the "Reply to" feature in the MyoVision software as Gmail will block emails sent



Note: If you need help at any point, please reach out to our technical support team at 1-800-969-6961, option 2. (International: 206-448-3464) Office hours are 8:00am - 4:00pm M-F, PST.

