

How Do I Disable Fast Startup?


Last Modified on 11/14/2018 9:37 am PST

Summary


Windows has a feature called "**Fast Startup**" that allows for your computer to boot faster. As a result, when the computer is shut down it does not register a restart. MyoVision Technical Support recommends turning this setting off.

To tun this setting off:


1. Open the Control Panel and search for "power button options".



2. Click "Change settings not currently available."



Left click the box next to "**Fast Startup**" to remove the check mark.



Click "Save Changes."

Now when you shut down your computer it will register it as a "Restart."

