

# Is MyoVision compatible with Windows 10?

Last Modified on 01/06/2020 3:46 pm PST

Short answer: Yes!

MyoVision WireFree System is fully compatible with Windows 10. If you have an M8000 system Windows is not fully compatible after Windows XP. Please see our support article for why. [End of Life for M8000 Wired System](#)

If you are buying a new system or upgrading your old one to Windows 10 see our support article which guides you through what you might need to know and how to install MyoVision. It includes tips on settings you might need as well as links to how to move your patient database and setting up your MyoVision software on Windows 10. [How to move your MyoVision system to a new computer](#)

---